

Non-invasive detection of Helicobacter pylori

HELIC® ABT Reader


Diagnostics made easy

- Real-time point-of-care testing
- Sensitivity 95 %, specificity 94,4 %
- Innovative breath test with the indicator tubes

HELIC ABT Reader allows for an automated air sampling and reading the color of the indicator tube. The indicator tube serves as the single-use sensitive element of the system, thus the Reader does not require periodical maintenance.



Thanks to the QR-code the software of the **HELIC ABT Reader** receives the batch parameters and then guides the operator throughout the testing procedure step-by-step to avoid human factor errors.



*Computer is not included

